


AeroEclipse^{*} XL

BAN^{*} Nebulizer

TRUETT MEDICAL
INTERNATIONAL


140-CA-EN

- EN** Sign up to receive a reminder when it's time to replace your nebulizer.
- DA** Tilmeld dig for at modtage en påmindelse, når det er tid til at udskifte forstøveren.
- DE** Melden Sie sich an und erhalten Sie eine Erinnerung, wenn Ihr Vernebler ausgetauscht werden sollte.
- EL** Εγγραφείτε για να λαμβάνετε μια υπενθύμιση όταν έρχεται η ώρα να αντικαταστήσετε τον νεφελοποιητή.
- ES (ES)** Insíbase para recibir un recordatorio cuando llegue el momento de sustituir el nebulizador.
- ES (LA)** Insíbase para recibir un recordatorio cuando llegue el momento de sustituir el nebulizador.
- ET** Registreeruge, et saada meeldetuletus, kui on aeg nebulisaator välja vahetada.
- FI** Tilaa muistutus nebulisaattorin vaihtamisesta.
- FR** Inscrivez-vous pour recevoir un rappel quand il sera temps de remplacer votre nébuliseur.
- IT** Registrati per ricevere un promemoria quando è il momento di sostituire il nebulizzatore.
- KO** 네불라이저 교체 시기에 알림을 받으려면 등록하십시오.
- LT** Užsiregistruokite, kad atėjus laikui pakeisti puršiklį gautumėte priminimą.
- LV** Reģistrējieties, lai saņemtu atgādinājumu, kad būs pienācis laiks nomainīt smidzinātāju.
- MS** Daftar untuk menerima peringatan apabila tiba masanya untuk menggantikan penebula anda.
- NL** Meld u aan om een herinnering te ontvangen wanneer het tijd is om uw vernevelaar te vervangen.
- NO** Registrer deg for å motta en påminnelse når det er på tide å skifte ut forstøveren.
- PL** Rejestracja umożliwiała otrzymywanie powiadomień o zbliżającym się terminie wymiany nebulizatora.
- PT** Inscreva-se para receber um lembrete quando chegar a altura de substituir o seu nebulizador.
- RO** Înscrieți-vă pentru a primi un memento atunci când este timpul să înlocuiți nebulizatorul.
- SI** ඔබේ නිහාරිකාරකය ආදේශ කිරීමට කාලය පැමිණි විට සිතිකැඳවීමක ලැබීමට ලියාපදිංචි වන්න.
- SL** Prijavite se za prejem opomnika o zamenjavi nebulizatorja.
- SV** Registrera dig för att få en påminnelse när det är dags att byta ut nebulisatorn.
- TA** உங்கள் நெபுலைசரை மாற்ற வேண்டிய நேரம் வரும்போது நினைவுட்டலைப் பெற பதிவு செய்க.
- TH** สมัครรับการแจ้งเตือนให้เปลี่ยนเครื่องพ่นละอองยาเมื่อถึงเวลา
- TR** Nebulizatörü değiştirme zamanı geldiğinde hatırlatma almak için kaydolun.
- VI** Đăng ký để nhận lời nhắc khi đến thời điểm cần thay máy xông khí dung.
- ZH (HK)** 請即登記接收通知，以便提醒您在適當時候更換您的霧化器。
- ZH (TW)** 註冊以接收何時須更換噴霧器的提醒。
- HE** הרשם כדי לקבל תזכורת כאשר הגיע הזמן להחליף את הנבולאיזר שלך.

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
1

INTENDED USE

The nebulizer is intended to be used to deliver aerosol drug to the patient's lungs. It is to be used by infant, child and adult patients. The patient is to be awake and can breathe unassisted. The environments of use include hospitals, clinics and the home. The nebulizer is single patient, reusable when used in the home. It is multi-patient, reusable when cleaned according to the hospital instructions.

Breath actuation only produces aerosol when the patient inhales. Aerosol is not produced when the patient exhales or takes a break in treatment. The clinical benefit is more drug is available to the lungs. Less drug is lost to the environment.

⚠ Caution: Breath actuation delivers drug on inspiration only. This may result in an increase of the dosage to the lungs. Health care providers should take this into consideration with regards to dosage.

If you have concerns regarding the use of the nebulizer, please speak with your health care provider.

Notes

- The nebulizer must only be used by a person who understands the instructions for use and can use the device safely. The treatment must be carried out by a responsible person if the patient is not able to use the nebulizer on their own. Infants, children, and persons with limited capabilities must be supervised.
- The nebulizer is designed to aerosolize solutions and suspensions.
- Only use with drugs prescribed by a health care provider.
- Do not use in an anaesthetic or a ventilator breathing system.

BEFORE EACH USE

Ensure these instructions have been read and are kept available at all times. Carefully examine the device. If there is any visible dust or other debris inside, then clean it. Check inside the **EZ Twist*** tubing (F) end for fraying or debris, clean it out if needed. Replace the nebulizer and the tubing (F) immediately if they are damaged or have missing parts.

⚠ Warning: Do not leave the nebulizer unattended with children. Contains small parts that can cause choking and tubing that can cause strangulation.

INSTRUCTIONS FOR USE

1 Ensure the exhalation valve (B) on the mouthpiece (A) faces down.

⚠ Warning: It is important to prevent exhaled breath from blowing near the eyes. Refer to the drug product monograph for more information.

① Note: For patients unable to use the mouthpiece, the nebulizer can be used with a **ComfortSeal*** Mask.

2 Inhale and exhale through the nebulizer mouthpiece (A) to verify the movement of the green feedback button (C) on the top of the nebulizer. The button should move up and down freely.

Notes:

- To operate the nebulizer in continuous mode (I), rotate the mode selector (G) clockwise to the continuous mode position .
- To return the nebulizer to breath actuated mode (H), rotate the mode selector (G) counterclockwise to the breath actuated position .

3 Unscrew and remove the top of the nebulizer. Place prescribed drug into the nebulizer cup (E). Reattach the top and gently hand tighten.

① Note: Keep the nebulizer in an upright position so the drug does not spill.

4 Use With Compressor: Attach one end of the tubing (F) to the fitting in the bottom of the nebulizer and the other end to the compressor. Ensure both ends are securely attached. Turn on the compressor.

① Note: The nebulizer may be used with an **Ombra*** Compressor or an air source capable of delivering 103 to 170 kPa (15 to 25 psi) and 4 L/min when connected to the nebulizer.

Hospital Use With Wall Air: Attach one end of the tubing (F) to the fitting in the bottom of the nebulizer and the other end to the flowmeter or compressed air source. Ensure both ends are securely attached. Turn on the compressed air source. Set the flowmeter to 7 to 8 L/min with a compressed air source capable of delivering 345 kPa (50 psi).

5 Place the mouthpiece (A) in mouth, and inhale slowly and deeply. As you inhale, the green feedback button (C) will move into the down position. This indicates the nebulizer is producing aerosol. Exhale normally. As you exhale, the green feedback button (C) will move into the up position. This indicates the nebulizer is not producing aerosol.

Notes

- Do not place the lips over the exhalation valve (B) on the bottom of the mouthpiece (A). Doing so will prevent the valve from functioning properly.
- Do not cover the inhalation ports (J) on the nebulizer. Doing so will reduce volume of inhaled air.
- Ensure the nebulizer remains upright for proper nebulization.

6 Continue to inhale from the nebulizer until it starts to sputter, or as instructed by your health care provider.

7 Turn off the compressed air source.

① Note: Disconnect the tubing (F) from the compressed air supply after use.

HOME USE – CLEANING

⚠ Caution: Follow the cleaning instructions. Failure to do so may affect drug delivery.

Nebulizer: The nebulizer (excluding the tubing (F)) can be cleaned after each use.

1 The mode selector (G) should be in the breath actuated position .

2 Remove the tubing (F) from the bottom of the nebulizer.

3 Unscrew and remove the nebulizer top (D).

4 Shake out excess drug.

5 Use any of the following methods to clean the nebulizer:

A Hand Washing

i. Soak the three parts in a mild solution of liquid dish detergent and lukewarm clean water for about 5 minutes.

ii. Agitate gently.

iii. Rinse thoroughly with clean water. Shake out excess water.

B Dishwasher

i. Secure the three parts in a basket on the top rack.

ii. Wash in a normal wash cycle with dishwasher detergent and a rinse aid if available.

iii. Shake out excess water.

① Note: Dishwashing with overly dirty dishes is not recommended.

6 Briefly connect the tubing (**F**) to the bottom of the nebulizer cup (**E**) and the compressed air supply to clear the air pathway.

7 Allow the parts to air dry thoroughly before reassembling.

8 Reattach the nebulizer top (**D**), making sure to line up the white cylinder with the stem of the nebulizer cup (**E**). Gently hand tighten until a stop is detected.

9 Reattach the mouthpiece (**A**) to the nebulizer with the exhalation valve (**B**) facing down.

10 After the nebulizer is fully dry and assembled, store in a clean plastic bag or container.

EZ Twist® Tubing: To remove condensation on the inside of the tubing (**F**), connect one end of the tubing (**F**) to the compressed air supply. Allow the air to run without the nebulizer attached. This will force air through the tubing, drying the inside. To clean the exterior, wipe with a damp cloth. Allow to air dry thoroughly.

HOME USE – DISINFECTION

The nebulizer (excluding the tubing (**F**)) can be disinfected each day.

Follow the cleaning instructions. Do not dry or reassemble the nebulizer.

1 Use any of the following methods to disinfect the nebulizer.

A Boiling

i. Place the three parts in a pot with boiling water and boil for 5 to 15 minutes.

① Note: Place a small rack on the bottom of the pot to prevent parts from coming in contact with the hot bottom.

ii. Remove from the water without touching the water.

iii. Allow the parts to air dry thoroughly before reassembling.

B Electronic Steam Disinfectors

i. Place the three parts in an electronic steam disinfectors. Follow the manufacturer's instructions for use.

ii. Remove the parts from the unit. Allow the parts to air dry thoroughly before reassembling.

C Microwave Steam Bag

i. Place the three parts in a microwave steam bag (such as Quick Clean® Micro-Steam® bag). Follow the manufacturer's instructions for use.

ii. Remove from the bag. Allow the parts to air dry thoroughly before reassembling.

D Isopropyl Alcohol

i. Soak the three parts in 70 % isopropyl alcohol for 5 minutes.

ii. Rinse the parts thoroughly with sterile water.

iii. Allow the parts to air dry thoroughly before reassembling.

E Bleach

i. Soak the three parts in a 1:50 bleach to water solution (1 tablespoon (15 mL) bleach in 3/4 cups (800 mL) water) for 3 minutes.

ii. Rinse the parts thoroughly with sterile water.

iii. Allow the parts to air dry thoroughly before reassembling.

① Note: You may substitute a respiratory equipment disinfectant (such as Milton®, Dodie®, Control III® or Cidex®) for the bleach solution. Follow the disinfectant instructions for use.

F Vinegar

i. Soak the three parts in a 1:3 vinegar to water solution (1 cup (250 mL) vinegar in 3 cups (750 mL) water) for 1 hour.

ii. Rinse the parts thoroughly with sterile water.

iii. Allow the parts to air dry thoroughly before reassembling.

2 After the nebulizer is fully dry and assembled, store in a clean plastic bag or container.

① Note: The nebulizer (excluding tubing) can be cleaned and disinfected up to 728 times.

HOSPITAL INSTRUCTIONS FOR MULTI-PATIENT RE-USE

⚠ Warning: Clean, disinfect and sterilize the nebulizer before using it with a different patient. This will prevent cross infection.

① Notes:

- The nebulizer (excluding the tubing (**F**)) can be cleaned, disinfected and sterilized up to 150 times.
- The tubing must be replaced between patients.

CLEANING AND DISINFECTION

1 The mode selector (**G**) should be in the breath actuated position .

2 Remove the tubing (**F**) from the bottom of the nebulizer.

3 Unscrew and remove the nebulizer top (**D**).

4 Shake out excess drug.

5 Soak the three parts in a 0.5 % solution of neodisher® MediClean forte (or similar) cleaning agent for 5 minutes, agitating gently.

6 Load the parts into a washer-disinfector compliant to ISO 15883-1. Ensure the parts are oriented upright.

7 Run the cleaning and disinfection cycle.

① Notes:

- Use cleaning agent neodisher® MediClean.
- The use of a rinse agent is optional.

8 Following disinfection, remove the parts from the washer-disinfector. Shake out excess water. Allow the parts to air dry thoroughly.

STERILIZATION

1 After cleaning, disinfection, and air dry, pack and seal the nebulizer parts in a sterilization pouch compliant to ISO 11607-1 & 2 or EN 868-5.

① Notes:

- The pouch should be no smaller than 70 mm x 120 mm.
- The white cylinder portion of the nebulizer top (**D**) should be positioned in the pouch so it is not bent or pressed into corners where it can become bent (Figure ).

2 Load the pouch in the sterilizer.

3 Sterilize according to ISO 17665-1 Moist Heat with any of the following parameter sets:

A

- 3 pre-vacuum phases
- 134 °C sterilization temperature
- 3 minutes hold time
- 1 minute dry time

B

- 3 pre-vacuum phases
- 121 °C sterilization temperature
- 20 minutes hold time
- 1 minute dry time

① Note: Do not exceed a temperature of 137 °C.

4 Reassemble the nebulizer prior to next use.

5 Store in a clean plastic bag or container.

① Notes

- No modifications of the nebulizer are permitted.
- The nebulizer parts comply with the relevant requirements of ISO 27427.
- The materials used for the components may not be compatible with emulsions. In such cases, information should be sought from the drug supplier.
- Not made or manufactured with bisphenol A (BPA), di(2-ethylhexyl) phthalate (DEHP), latex or lead.
- This nebulizer can be disposed of with domestic waste unless this is prohibited by the disposal regulations prevailing in the respective member countries.
- Any serious incident that has occurred in relation to the device should be reported to the manufacturer and the competent authority of the Member State, the MHRA or the Saudi Food and Drug Authority, as appropriate.

Type	Pneumatically powered jet nebulizer
Breath Actuation Operation	Triggered mechanically (by the patient on inhalation) by an inhalation flow over 12 L/min
Driving Gas	Air
Driving Gas Pressure – Compressor (at nominal 4 L/min flow when connected to the nebulizer)	Maximum: 172 kPa (25 psi) Nominal: 140 kPa (20 psi) Minimum: 103 kPa (15 psi)
Driving Gas Flow – Wall Air (at nominal 345 kPa (50 psi) pressure when connected to the nebulizer)	Maximum/Nominal: 8 L/min Minimum: 7 L/min
Maximum Pressure for Tubing at 22 °C ± 2 °C Ambient Conditions	345 kPa (50 psi)
Maximum Fill Volume	6 mL
Maximum A-Weighted Sound Pressure Level	63 dBA
Mass of Nebulizer (without tubing)	0.03 kg (0.07 lb)
Operating Temperature Limits	+15 °C (59 °F) to +40 °C (104 °F)
Operating Relative Humidity Range	15 % to 95 % RH
Storage Temperature Range	-40 °C (-40 °F) to +60 °C (140 °F)
Storage Relative Humidity Range	Ambient to 95 % RH
Recommended Replacement Period	6 months

COMPRESSOR — ISO 27427 PERFORMANCE DATA	At 3.5 L/min (100 kPa)	At 4 L/min (130 kPa)	At 4.5 L/min (160 kPa)
Aerosol Output ¹	10.1 mg	10.8 mg	11.5 mg
Aerosol Output Rate ¹	1.7 mg/min	1.9 mg/min	2.2 mg/min
Percentage of Fill Volume Emitted Per Minute ¹	3.5 %	3.9 %	4.4 %
Residual Volume ¹	1.2 mL	1.1 mL	1 mL
Mass Median Aerodynamic Diameter (MMAD)	4.0 µm ²		

COMPRESSOR — ISO 27427 PERFORMANCE DATA	At 3.5 L/min (100 kPa)	At 4 L/min (130 kPa)	At 4.5 L/min (160 kPa)
Geometric Standard Deviation (GSD)	2.1 ²		
Particle size % > 5 µm	39 % ²		
Particle size % 2 µm to 5 µm	41 % ²		
Particle size % < 2 µm	20 % ²		
Respirable Fraction (5 µm and smaller)	61 % ²		

WALL AIR — ISO 27427 PERFORMANCE DATA	At 7 L/min	At 8 L/min
Aerosol Output ¹	14.8 mg	16.1 mg
Aerosol Output Rate ¹	3.4 mg/min	3.9 mg/min
Percentage of Fill Volume Emitted Per Minute ¹	6.9 %	7.8 %
Residual Volume ¹	1 mL	1 mL
Mass Median Aerodynamic Diameter (MMAD)	3.4 µm ³	
Geometric Standard Deviation (GSD)	2.3 ³	
Particle size % > 5 µm	33 % ³	
Particle size % 2 µm to 5 µm	40 % ³	
Particle size % < 2 µm	27 % ³	
Respirable Fraction (5 µm and smaller)	67 % ³	

¹ 2 mL fill volume of sodium fluoride 2.5 % (M/V) solution. TMI data on file.

² 2.5 mL fill volume of salbutamol 0.1 % (M/V) solution. TMI data on file.

³ 3 mL fill volume of sodium fluoride 2.5 % (M/V) solution. TMI data on file.

① Notes

- A suspension and/or high viscosity solution can alter the particle size distribution curve, the mass median aerodynamic diameter (MMAD), aerosol output and/or aerosol output rate which then can be different from those disclosed in these instructions.
- Disclosures for nebulizer performance are based upon testing that utilizes adult ventilatory patterns and are likely to be different from those stated for paediatric or infant populations.

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