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Aerobika^{*} ARX^{*}

Connected Technology

TRUETT MEDICAL
INTERNATIONAL



FIGURE 1A

***Aerobika*^{*} OPEP Device**
Dispositif de PEPO *Aerobika*^{*}




FIGURE 1B

Arx* Connected Technology (Arx* Device)
Technologie Connectée Arx* (Appareil Arx*)

Arx* Electronic Module (electronic module)
Module Electronique Arx* (module électronique)

Adaptor / Adaptateur

EN

- A Module port
- B Arx* device label and serial number
- C Button
- D Light ring
- E USB-C port
- F Adaptor port
- G Membrane
- H Charging cord

FR

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Membrane Loosened /
Membrane desserrée

INTENDED USE

The **Aerobika® Arx®** device is an accessory to the **Aerobika®** OPEP device and is designed to track and provide feedback on **Aerobika®** OPEP device usage.

HOW THE ARX® DEVICE WORKS

- The device provides real-time feedback during **Aerobika®** OPEP usage.
- The device collects information related to airway clearance therapy during usage.
- The information is transmitted via Bluetooth® to your smartphone.
- Use the **AeroLiv®** mobile app to view information collected by the device.

HARDWARE AND SOFTWARE REQUIREMENTS

The **Arx®** device can be operated with the **AeroLiv®** mobile app on both iOS and Android operating systems. The applications are available in Apple App store and Google Play store. The applications work on version iOS 16 and Android 10 or above. Communication of the **Arx®** device to the **AeroLiv®** Application is via Bluetooth® 5.0 (BT LE) and requires a mobile device that accommodates this requirement. We do not recommend the use of applications on tablets and iPads.

For additional information on the use and benefits of the **AeroLiv®** mobile app visit aeroliv.ca.

INDICATION(S) AND CONTRAINDICATIONS

For indications and contraindications for the **Aerobika®** OPEP device, consult the Instructions for Use.

Single patient use only.

USE INSTRUCTIONS

Before use, ensure these instructions, and those for the **Aerobika®** OPEP device, have been read and are kept available at all times. If you have questions about the performance or use of this product, please contact your healthcare professional prior to use. For device part descriptions see Figure 1.

GETTING STARTED

⚠ Warning: Do not leave the device unattended with children. Contains small parts that can cause choking.

Part 1: Assembly

1 Remove all the pieces from the package. See Figure 1 for a full list.

2 Remove the mouthpiece from the **Aerobika®** device.

3 Attach the mouthpiece to the Adaptor. It will only fit on one end.

4 Attach the **Aerobika®** device to the other end of the Adaptor. See Figure 2 to confirm you have assembled your device correctly.

ⓘ Note: The device can be used directly out of the package.

Part 2: Activate Module

1 To turn on for the first time, press and hold the button on the top of the Module until the light ring lights up. This will take approx. 3 seconds.

2 Let go of button once lights appear.

- If the device doesn't light up after 10 seconds, the battery may be dead. You will need to charge your device before you can record your first airway clearance therapy session.

3 The device is now "awake" and is ready to use when the front of the ring pulses a soft blue light.

- If the device displays an amber pulse at the front of the light ring, the device's battery needs to be charged. You can still perform an airway clearance therapy session, and the device will track your performance. However, to conserve the battery, it will not provide visual feedback.

4 The device will go to sleep after idling for 30 seconds and the blue light will turn off. When the device is picked up or moved, the blue light will pulse to indicate it is "awake" again.

5 Your **Arx®** device is now ready to use. It will record several days of therapy data but will need to be paired with your phone to download data. To pair with your phone, please see the next section.

ⓘ Note: The **Aerobika®** device will still work even if the Electronic Module battery is dead, however the treatment will not be recorded.

Part 3: Using your Arx® Connected Technology with the AeroLiv® Mobile App

AeroLiv® is a mobile application available on both Apple and Android devices. To learn more about the app and to begin using it with your **Arx®** device visit aeroliv.ca or scan the QR code below with your phone camera.

COMPLETING AN AIRWAY CLEARANCE THERAPY SESSION

Prior to use check the following:

- Carefully examine the device. If there is any visible dust or other debris inside, then clean the device (see Cleaning section, pg. 5).
- Remove the Electronic Module from the Adaptor and check the membrane for tears, puncture holes or other damage. Replace immediately if it is damaged or if device has missing parts.

Aerobika® Arx® Therapy Settings

The **Arx®** device includes therapy settings that can be changed in the **AeroLiv®** app.

⚠ Caution: Prior to use check that all therapy settings are correct. Changes to settings must be determined by a healthcare professional.

- The device has a default treatment time of 10 minutes, but this can be changed in the **AeroLiv®** app.
- If you track your device-use by counting your breaths in reps and sets, this can also be programmed in the **AeroLiv®** app.
- Always use the method prescribed by your healthcare professional.

① Note: Changing a therapy setting while in an airway clearance therapy session will start a new session and restart your treatment time.

Using the Aerobika® Arx® Device

1 The assembled device can be used in any orientation. It is recommended that treatments be done in a comfortably seated position.

2 Pick up the **Aerobika® Arx®** device. The device will automatically calibrate each time it “wakes up” so wait for the blue pulse before beginning your treatment.

- If the blue pulse does not appear, your device may be in Travel Mode, or the battery may require charging. For more information about travel mode see the Travel Mode section, pg 6.
- If you start to exhale into the device before the light pulse reappears, your device has not had time to calibrate and will not record your airway clearance therapy session. Place the device back on a flat surface and wait for the light pulse to reappear, this should only take 1-2 seconds.

3 Follow the steps included in the **Aerobika®** device instructions for use for additional guidance on completing an airway clearance treatment.

4 The goal of each breath is to see green on the light ring. This indicates you are exhaling through the device within the desired range.

- If you are exhaling too softly, a solid blue light will appear in the front of the light ring. To correct your technique, exhale more forcefully with a firm steady breath until the light ring turns green.
- If you are exhaling too hard, the light ring will flash red. To correct your technique, exhale softer until the light ring turns green.

5 When you see the green lights move in a circular pattern you have reached your treatment time or breath-count goal. The pattern will repeat after each additional exhalation.

Airway Clearance Therapy Technique Guide

Solid Blue	You are exhaling too softly. Exhale harder.
Solid Green	Any combination of green lights means you are in the target zone. This is your goal for each breath.
Flashing Red	You are exhaling too hard. Exhale softer.
Circular Green	You have reached your treatment time or breath-count goal

See the Device Light Pattern Guide for more detail.

① Note: The **Arx®** device assumes that any exhalation that occurs within a 30-minute period is part of the same treatment session. If you wait longer than 30 minutes between exhalations, the device will record it as the start of a new session.

① Note: Changes in altitude may impact device function. Remove and reinstall the Electronics Module from the Adaptor while at altitude to ensure proper function.

Using your Aerobika® OPEP Device with Arx® Device with your Nebulizer

The device may also be used with a small volume nebulizer with a 22 mm fitting. Your healthcare professional will advise which medication to use for combined treatments.

Attach the nebulizer to the back of the **Aerobika®** device as described in the **Aerobika®** device Instructions for Use. Follow the instructions provided with your **Aerobika®** device for how to use with your nebulizer.

CLEANING

- The Electronic Module contains electronic components and CANNOT be put in water, boiled, microwaved, or steamed.
- Proper and regular cleaning of the **Arx®** device is important for correct operation.

A Disassembly

You must separate the **Aerobika®** device and **Arx®** device into components before cleaning.

- Remove the **Aerobika®** body from the Adaptor.
- Hold the Adaptor, with the Mouthpiece facing away from you.
- To remove the Electronic Module, grip firmly and rotate counterclockwise while pulling away from the Adaptor (see Figure 3).
- Remove the Mouthpiece from Adaptor.
- The green silicone membrane should be loosened from its attachment point but remain attached to the Adaptor during cleaning (see Figure 4).
- Follow the disassembly instructions included with your **Aerobika®** device.

B Cleaning Procedure

i. The Adaptor should be cleaned whenever you clean your **Aerobika*** device. It can be cleaned by handwashing in a warm water/liquid dish detergent solution, rinsed with clean water and shaken to remove any excess.

ii. The outer surface of the Electronic Module can be wiped with a damp (not wet) alcohol swab and air dried.

- To create an alcohol swab, add 70% isopropyl alcohol to a sterile cotton swab or use a sealed store-bought alcohol wet wipe.

iii. Ensure the Adaptor is completely dry before reattaching the Electronic Module.

iv. Follow the cleaning instructions included with your **Aerobika*** device.

C Complete Assembly of Arx* and Aerobika* devices after cleaning

i. Reattach the Mouthpiece to the front of the Adaptor.

ii. Attach the **Aerobika*** body to the other end of the Adaptor.

iii. Reseat the green silicone membrane you loosened prior to cleaning on the circular port on the top of the Adaptor.

iv. Reattach the Electronic Module following the steps in Figure 3 in reverse order. Align the Electronic Module to the Adaptor port and attach it to the Adaptor. The two pieces will not fit together unless you have the Electronic Module angled correctly.

v. Once the two circular ports are engaged, gently rotate the Electronic Module clockwise until it is in-line with the Adaptor, and it is firmly locked in place. Do not overtighten.

vi. The device is now ready to use.

CHARGING YOUR DEVICE

If your device needs to be charged, the light ring will turn amber and show the Low-Power pattern.

1 Plug any end of the USB C charging cord into the port on the back of the Module.

2 Plug the other end of the cord into a power source (USB C wall adaptor not included).

 **Warning:** Use only with an IEC 60601-1 certified power source.

The front of the light ring will pulse green twice to indicate the device is charging, then default to the blue Standby pulse while it's awake.

When the device falls asleep, it will pulse green. The front of the light ring will change to a solid green once charging is complete.

Do not charge while in use.

OPTIONAL WIRELESS CHARGING (CHARGING PAD NOT INCLUDED)

With the Adaptor still attached, place the Module button-side down in the center of your charging pad. The front of the light ring will display the charging pattern when it is placed in the charging area.

If the ring does not light up, reposition on the charging pad until the light is visible.

TRAVEL MODE

When you move or handle your **Arx*** device, it “wakes up” and gets ready for use. Travel mode powers down the device and disables Bluetooth® to prevent the device battery from draining when it is being moved. The device won't record a therapy session or display any lights when in travel mode.

 **Note:** Use the **AeroLiv*** app to place your **Arx*** device into travel mode.

To exit travel mode, press the button on top of the Module with a single click. When the Module exits travel mode, the light ring will display a Wake-Up pattern and will then pulse blue when it is ready to use.

TROUBLESHOOTING

Please contact support at aeroliv@trudellmed.com if any of the following occur:

- Your device does not turn on.
- The light does not light up when plugged-in to charge.
- Your device continues to display the error pattern after a hard reset.

How To Perform a Hard Reset

Although rare, you may need to perform a reboot of your Module. If you are not using the device, and the light ring is pulsing red for more than 30 seconds, a hard reset is required.

1 Press and hold the button on the **Arx*** device down until the Wake-Up pattern appears (10-15 seconds).

 **Note:** Ensure you hold the button firmly. You should only hear the button click once when you press down on it, and again when you let go after seeing the Wake-up pattern. If you have held the button for more than 20 seconds without seeing the Wake-Up pattern, stop pressing the button and try again.

2 When the device has rebooted the device will display a blue pulse when it is ready to use.

 **Caution:** All therapy settings will be set to their default values after a hard reset.

AEROBIKA* ARX* DEVICE LIGHT PATTERN GUIDE

Colour	Pattern Name	Pattern	Meaning
Blue	Ready-to-Go	The front of the light ring displays a soft blue pulse.	The device is ready to use. It should always light up blue prior to starting your airway clearance therapy.
	Breathe-Harder	The front of the light ring displays a steady blue while you are blowing into the device.	This pattern appears during airway clearance therapy. If you see this, your exhalation pressure is too low, and you need to blow harder through the device until the colour turns green.
	Bluetooth®-Pairing	Blue lights sweep up the light ring and then back down, multiple times in a row.	When device pairing is initiated by the AeroLiv* app, the Arx* device will respond by showing this pattern. While the pattern appears, the device is attempting to set up a Bluetooth® connection. The pattern will repeat until a successful connection is made or until 30 sec has elapsed (whichever comes first).
Green	Aerobika*	The lights on the light ring display green while you are blowing into the device. The lights will light up from the front of the device to the back with increasing exhale pressure.	This pattern appears during airway clearance therapy. If you see this, your exhalation pressure is within the desired range. Good job!
	Therapy-Complete	Green lights circle the light ring in a clockwise direction, repeating three times total.	You have reached your targeted airway clearance therapy treatment time or breath-count goal. The pattern will repeat after each additional exhalation.
	Charging	The front of the light ring displays a green pulse until device is fully charged.	Device is charging but not yet charged.
	Charging-Complete	The front of the light ring displays a steady green until the device is unplugged.	The device is now fully charged.
Red	Slow-Down	The entire light ring flashes red two times per second and continues while the pressure is too high.	This pattern appears during airway clearance therapy. If you see this, your exhalation pressure is too high exhale less forcefully until the pattern changes and the colour turns green.
	Error	The entire light ring flashes red three times followed by a two second pause.	The device has encountered an error and is trying to fix it. If the ring light continues to flash for more than 30 seconds, it may require a hard reset. See "Hard Reset" instructions.
Multi-Coloured (Pink to Purple to Blue)	Wake-Up	Multi-coloured lights sweep up the light ring and then back down, three times in a row.	This pattern indicates the device is powering on and occurs when turning on for the first time, when exiting travel mode, and after a hard rest.
	AeroLiv* -Connection	Multi-coloured lights sweep up the light ring from front to back, repeating three times.	If you were attempting to pair your device to the app, this pattern indicates that the Bluetooth® pairing was successful.
	Firmware Updated	Multi-coloured lights sweep up the light ring from front to back, repeating three times.	If you accepted a firmware update, this pattern indicates that the update is complete, and your device is ready to use. This pattern will not appear if your firmware update was scheduled for a later time.
Amber / Orange	Low-Power Pattern	The front of the light ring displays an amber pulse.	The device has low battery power and needs to charge. You can still perform a therapy session, and the device will track your performance but will not provide visual feedback to conserve the battery.

⚠ CAUTIONS

Do Not Wet Electronics: Do not get the Module wet or submerge it in liquid as it can malfunction due to liquid damage.

Dizzy or Lightheaded: If you feel dizzy or lightheaded, or have any other discomfort during treatment, discontinue therapy and contact your physician or healthcare provider.

⚠ WARNINGS

Small Parts: Do not leave the device unattended with children. It contains small parts that can cause choking.

Proximity to RF Peripherals: Portable RF communications equipment (including peripherals such as antenna cables and external antennas) should be used no closer than 30 cm (12 inches) to any part of the Module, including cables specified by the manufacturer. Otherwise, degradation of the performance of this equipment could result.

Use of accessories: Use of accessories, transducers, and cables other than those specified or provided by the manufacturer of this equipment could result in increased electromagnetic emissions or decreased electromagnetic immunity of this equipment and result in improper operation.

Use adjacent to or stack with other equipment: Use of this equipment adjacent to or stacked with other equipment should be avoided because it could result in improper operation. If such use is necessary, this equipment and the other equipment should be observed to verify that they are operating normally.

Do Not Open: Opening the **Arx*** device by unscrewing the screws will void the device warranty and should not be attempted. If you require assistance with your device, contact customer service at: trudellmed.com/ca/en-CA/contact-us

No modification of this equipment is allowed.

STORAGE AND DEVICE INFORMATION

After the device is fully dry and assembled, store in a clean plastic bag or container.

If a replacement Adaptor, membrane, or charging cable is required, it should be purchased from www.aeroliv.ca.

The intended environment of use is Home Healthcare Environment. The patient is an intended operator of the device.

The **Aerobika*** OPEP device with **Arx*** connected technology needs to accurately measure pressure between 0-50 cm H₂O. The device has been qualified to remain functional against the presence of electromagnetic disturbances up to the values specified in table 2, but exposure above those indicated could lead to an error in pressure measurement. Do not use in the presence of devices generating high electromagnetic fields such as magnetic resonance imaging (MRI) equipment.

The **Arx*** device consists of an integrated radio module with the following characteristics:

- BLE radio module
- Frequency band: 2.4-2.4835 GHz
- Modulation: GFSK
- E.I.R.P.: 2.8 dBm
- Occupied bandwidth: 1.05 MHz

DISPOSAL

The **Arx*** device should be replaced three years after first use. The adaptor and membrane should be replaced every year.

Dispose of the **Arx*** device in accordance with the national regulations in force in your country.

Do not place used devices or batteries in normal waste containers. Contact your recycling company for this purpose. The device may be handed over to the manufacturer, distributor, or recycler. If the battery inside the device is damaged, the product must be returned to the manufacturer.

The **AeroLiv*** app can be used to erase the data on the device prior to disposal.

ADDITIONAL INFORMATION

Use Temperature Range	+15°C to +40°C
Storage Temperature Range	-20°C to +55°C
Use Humidity Range	30% to 70% Relative Humidity
Storage Humidity Range	10% to 90% Relative Humidity
Atmospheric Pressure Range	70 kPa to 106 kPa
Measurement Error	0 -50 cm H ₂ O = +/- 5% or +/- 0.5 cm H ₂ O, whichever is greater
Typical Operating Use Time	20 Minutes
Typical Service Life of Arx* Device	3 Years
Typical Service Life of Adapter and Membrane	1 Year
Typical Service Life of Aerobika* OPEP Device	Consult the Instructions for Use

- Product is free of Bisphenol

INFORMATION ON COMPLIANCE

The **Arx*** device was tested against all applicable tests in IEC 60601-1-2:2014+A1:2020 and found in compliance with all tests.

All applied tests are listed in the tables below.

Table 1: Emissions

Emission tests	Compliance
Radiated RF emissions CISPR 11	Group 1, class B
Conducted RF emissions CISPR 11	Group 1, class B
Harmonic emissions IEC 61000-3-2:2014	Class A
Voltage fluctuations and flickers emissions IEC 61000-3-3:2013	Complies

Table 2: Immunity

Immunity tests	Standard test level	Compliance
Radiated RF EM field IEC 61000-4-3:2006 + A1:2007 + A2:2010	10 V/m (80 – 2700 MHz), 80 % AM modulated (1 kHz)	10 V/m (80 – 2700 MHz), 80 % AM modulated (1 kHz)
Proximity fields from RF wireless communications equipment IEC 61000-4-3:2006 + A1:2007 + A2:2010	See table 3 below	See table 3 below
Conducted disturbances induced by RF fields IEC 61000-4-6:2013	3 VRMS (unmodulated) 0.15–80 MHz, 6 VRMS (unmodulated) in ISM and amateur radio bands between 0.15–80 MHz 80 % AM (1 kHz)	3 VRMS (unmodulated) 0.15–80 MHz, 6 VRMS (unmodulated) in ISM and amateur radio bands between 0.15–80 MHz 80 % AM (1 kHz)
Electrostatic discharge IEC 61000-4-2:2008	± 8 kV (contact) ± 2, 4, 8, 15 kV (air)	± 8 kV (contact) ± 2, 4, 8, 15 kV (air)
Surge IEC 61000-4-5:2014 + A1:2017	± 1 kV (line to line) 1.2/50 (8/20) Tr/Th ms	± 1 kV (line to line) 1.2/50 (8/20) Tr/Th ms
Electrical fast transients/bursts IEC 61000-4-4:2012	± 2 kV (peak), 5/50 Tr/Th ns, 100 kHz (repetition rate)	± 2 kV (peak), 5/50 Tr/Th ns, 100 kHz (repetition rate)
Voltage dips and interruptions IEC 61000-4-11:2004 + A1:2017	Voltage dip : 0 % for 0.5 cycles Voltage dip: 0% for 1 cycle Voltage dip: 70% for 25/30 cycles Short Interruption : 0 % during 250/300 cycles	Voltage dip : 0 % for 0.5 cycles Voltage dip: 0% for 1 cycle Voltage dip: 70% for 25/30 cycles Short Interruption : 0 % during 250/300 cycles
Power-frequency magnetic field IEC 61000-4-8:2009	30 A/m (50 Hz)	30 A/m (50 Hz)
Proximity magnetic fields IEC 61000-4-39:2017	See table 4 below	See table 4 below

Table 3: Proximity fields from RF wireless communications equipment specification

Test frequency (MHz)	Band (MHz)	Service	Modulation	Immunity Test Level ² (V/m)
385	380 – 390	TETRA 400	Pulse modulation, 18 Hz	27
450		GMRS 460, FRS 460	FM, ± 5 kHz deviation, 1 kHz sine ¹	28
710				
745	704 – 787	LTE Band 13, 17	Pulse modulation, 217 Hz	9
780				
810				
870	800 – 960	GSM 800/900, TETRA 800, iDEN 820, CDMA 850, LTE Band 5	Pulse modulation, 18 Hz	28
930				
1720				
1845	1700 – 1990	GSM 1800, CDMA 1900, GSM 1900, DECT, LTE Band 1, 3, 4, 25, UMTS	Pulse modulation, 217 Hz	28
1970				
2450				
5240	2400 – 2570	Bluetooth®, WLAN, 802.11 b/g/n, RFID 2450, LTE Band 7	Pulse modulation, 217 Hz	28
5500				
5785				
	5100 – 5800	WLAN 802.11 a/n	Pulse modulation, 217 Hz	9

Table 4: Proximity magnetic field specifications

Test Frequency	Modulation	Immunity Test Level (A/m)
30 kHz ¹	CW	8
134.2 kHz	Pulse Modulation ² 2.1 kHz	65 ³
13.56 MHz	Pulse Modulation ² 50 kHz	7.5 ³

¹ This test is applicable only to ME EQUIPMENT and ME SYSTEMS intended for use in the HOME HEALTHCARE ENVIRONMENT

² The carrier shall be modulated using a 50% duty cycle square wave signal.

³ r.m.s., before modulation is applied.

ISED NON-INTERFERENCE DISCLAIMER

This device contains license-exempt transmitter(s)/ receiver(s) that comply with Innovation, Science and Economic Development Canada's license-exempt RSS(s). Operation is subject to the following two conditions:

- 1** This device may not cause interference.
- 2** This device must accept any interference, including interference that may cause undesired operation of the device.

This device complies with the Canadian ICES-003 Class B specifications. CAN ICES-003(B) / NMB-003 (B).

FCC COMPLIANCE STATEMENT

This equipment has been tested and found to comply with the limits for a Class B digital device, pursuant to part 15 of the FCC Rules. These limits are designed to provide reasonable protection against harmful interference in a residential installation. This equipment generates, uses and can radiate radio frequency energy and, if not installed and used in accordance with the instructions, may cause harmful interference to radio communications. However, there is no guarantee that interference will not occur in a particular installation. If this equipment does cause harmful interference to radio or television reception, which can be determined by turning the equipment off and on, the user is encouraged to try to correct the interference by one or more of the following measures:

- Reorient or relocate the receiving antenna.
- Increase the separation between the equipment and receiver.

- Connect the equipment into an outlet on a circuit different from that to which the receiver is connected.
- Consult the dealer or an experienced radio/TV technician for help.

This device complies with Part 15 of the FCC Rules. Operation is subject to the following two conditions:

- 1** this device may not cause harmful interference, and
- 2** this device must accept any interference received, including interference that may cause undesired operation. Please note that changes or modifications not expressly approved by the party responsible for compliance could void the user's authority to operate the equipment.

Changes or modifications not expressly approved by the party responsible for compliance could void the user's authority to operate the equipment.

ISED/FCC RF EXPOSURE (SAR) STATEMENT

This equipment that is intended to be operated close to the human body is tested for body-worn Specific Absorption Rate (SAR) compliance. The SAR limit set by the FCC and ISED is 1.6 W/kg when averaged over 1 g of tissue. When carrying the product or using it while worn on your body, maintain a distance of 50 mm from the body to ensure compliance with RF exposure requirements. This equipment complies with ANSI/IEEE C95.1-1999 and are tested in accordance with the measurement methods and procedures specified in OET Bulletin 65 Supplement C.

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DECLARATION OF CONFORMITY TO HUMAN EXPOSURE LIMITS TO RADIOFREQUENCIES (DAS) OF ISDE/FCC

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MD Medical device
Dispositif médical

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 Importer
Importateur

 Manufacturer
Fabricant

 Date of manufacture
Date de fabrication

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 Information
Renseignements

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 Warning
Avertissement

 Refer to instruction manual
Consulter le manuel d'instructions

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 Designed and developed in Canada
Conçu et mis au point au Canada

MAX. 5V Device maximum voltage limit of 5 volts
Limite de tension maximale du dispositif de 5 volts

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Scan or contact Trudell Medical International to obtain a glossary of the symbols used on packaging and instructions for use.

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TRUDELL MEDICAL
INTERNATIONAL

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Trudell Medical International

725 Baransway Drive, London, Ontario,
Canada N5V 5G4

info@trudellmed.com www.trudellmed.com

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* trade marks and registered trade marks of Trudell Medical International (TMI). ® trade marks and registered trade marks of their respective companies.
© TMI 2025. All rights reserved. | * Marques de commerce et marques déposées de Trudell Medical International (TMI). ® Les marques de commerce