TRU ZONE*

PEAK FLOW METER

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ESPIRÓMETRO

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
479-GLO-EN

TRUDELL MEDICAL
INTERNATIONAL

People who use a valved holding chamber with their inhaler have **better control**.¹

Guidelines recommend chambers for children and adults.⁴

Ask your healthcare provider which chamber is right for you.

AeroChamber
Plus **Flow-Vu**[®]

Anti-Static Valved Holding Chamber

www.aerochamber.com

¹Levy ML, Hardwell A, McKnight E, Holmes J. Asthma patients' inability to use a pressurized metered-dose inhaler (pMDI) correctly correlates with poor asthma control as defined by the Global Initiative for Asthma (GINA) strategy: a retrospective analysis. *Prim Care Respir J*. 2013;Dec;22(4):406-11. ²Plaza V, Sanchis J. Medical personnel and patient skill in the use of metered dose inhalers: a multicentric study. *Respiration* 1998;65:195-8. ³Singh D et al. Effect of **AeroChamber Plus*** on the lung and systemic bioavailability of beclomethasone dipropionate/formoterol pMDI. *British Journal of Clinical Pharmacology* 2011;72(6):932-9. ⁴Global Initiative for Asthma. *Global strategy for asthma management and prevention 2017*. Available at: www.ginasthma.org/guidelines-gina-report-global-strategy-for-asthma.html.

EN	Instructions / Daily Record / ColorZone* Tapes.....	2
FR	Instructions / Enregistrement quotidien / Bandes ColorZone* .	4
ES	Instrucciones / Registro Diario / Cintas ColorZone*	7
VI	Hướng dẫn / Phiếu ghi nhận hàng ngày / Dải băng ColorZone*	9
AR	التعليمات / السجل اليومي / أشرطة ColorZone*	13

TRUZONE[®] PFM

Peak Flow Meter

Truly Easy

EN

⚠ Please read the instructions carefully before using.

Your **TruZone**[®] PFM has been designed for single patient use only.

INTRODUCTION


A peak flow meter is used to measure the fastest rate you can force air from your lungs. With asthma, your peak flow reading may drop before you even experience any symptoms. Tracking your readings will help detect these changes and may help you maintain control of your asthma.

TruZone[®] PFM can be used by either children or adults as part of an overall asthma management program, and must always be used as instructed. If you have questions about how to use this device or do not fully understand the user instructions, please contact your healthcare professional. **If you record an unusual reading, seek professional advice immediately.**

USING YOUR PEAK FLOW METER

It will be easiest for you to remember to record your peak flow reading if you associate it with something you do twice daily – like brushing your teeth or prior to taking your asthma medications. It is generally advised that peak flow be measured at 8:00 AM and 8:00 PM, but always as recommended by your healthcare professional. Your **TruZone**[®] PFM should be cleaned before first use and then weekly.

1 Visually inspect the device and make sure there are no defects or foreign objects lodged inside. This can be done by lightly blowing into the PFM to see if the yellow indicator moves, then resetting the indicator. Replace **IMMEDIATELY** if any defect is noticed.

2 To reset the yellow indicator, gently shake the device vertically until it is positioned directly behind the reset diamond.

3 Hold your **TruZone**[®] PFM by the clear body, ensuring the end of the device is not blocked.


4 Stand or sit up as straight as possible, while still being comfortable.

5 Take as deep a breath as possible. Do not inhale through the device.

6 Put the **TruZone**[®] PFM mouthpiece into your mouth – making sure to seal your lips tightly around the green mouthpiece. Don't block the mouthpiece with your teeth or tongue as this may result in an inaccurate reading.

7 Blow out as hard and as fast as possible.

The force of the air coming out of your lungs will cause the yellow indicator to move along the numbered scale.


① Note: Avoid coughing or spitting into the device to ensure accuracy.


8 Repeat steps 2 through 7 three times (or as instructed by your healthcare professional), resetting the yellow indicator between each reading.

9 Record the best of your three readings in your Daily Record.

USING YOUR PEAK FLOW METER WITH THE OPTIONAL COLORZONE[®] TAPE

ColorZone[®] Tapes are enclosed to help define your individual target zones. These are similar to the colours of a traffic light. Green signals that your asthma is under control and that your airways are open. Yellow signals that your condition may be worsening. You will likely start to experience symptoms, and the closer you move to the red zone the more pronounced the symptoms will become. The red zone indicates that you are in severe respiratory distress and immediate danger – **SEEK IMMEDIATE MEDICAL CARE!**

To apply, simply remove one piece of multicolored tape from the sheet by cutting along the dotted lines. To apply, peel the backing from the green zone. Place the tape in the middle of the scale and align the arrows below the word 'BEST' with your best peak flow number. Press the green zone into place. Cut the red zone just above the reset diamond (see illustration), peel off the remaining backing and press into place.

GREEN | Asthma Under Control

- 80-100 percent of personal best; asthma is well controlled

YELLOW | Worsening Condition

- High Yellow: 65-79 percent of personal best; asthma is getting worse
- Low Yellow: 50-64 percent of personal best

RED | Take Action

- Below 50 percent of personal best; **DANGER – SEEK IMMEDIATE MEDICAL HELP**

Once you have set the **ColorZone*** Tapes, they **SHOULD NOT** be moved unless you are advised to do so by a healthcare professional. Store the extra tapes in a safe place for use later.

HOW TO USE YOUR DAILY RECORD

Your Daily Record should be used to track your peak flow readings between healthcare appointments. There are spaces across the top to track up to 35 days, and an area below each date to track up to three separate peak flow readings. You have a choice of two methods when charting your peak expiratory flow results:

- **TruZone* PFM Method:** The scales on both your **TruZone*** PFM and Daily Record are identical, so it is easy to record your best peak flow number. Simply lay your **TruZone*** PFM directly onto the Daily Record. Align the device to the upper and lower limits of the scales of the Daily Record. Now mark an "X" next to the point where the yellow indicator rests, ensuring you are below the correct Date and Time on the graph.
- **Conventional Method:** Simply locate the point on the Daily Record, below the correct Date

and Time, corresponding to the best peak flow reading you wish to record. Mark an "X" at that point.

By tracking your peak flow reading twice each day, you can help your healthcare professional understand how well your asthma is controlled. Remember to take your Daily Record with you to every appointment.

Regardless of your peak flow reading, if you feel you are experiencing symptoms of an asthma attack seek immediate medical attention.

CLEANING YOUR TRUZONE* PFM

Your **TruZone*** PFM should be cleaned weekly.

1 Soak the device in warm soapy water, agitating gently. Do not use a brush or other object to clean the inside.

2 Rinse gently with clean water.

3 Gently shake out excess water and allow to air dry thoroughly in a vertical position.

4 You may also safely clean your **TruZone*** PFM using the normal cycle on the top rack of the dishwasher. Cleaning with overly dirty dishes is not recommended.

⚠ CAUTION: Never disassemble your **TruZone*** PFM or place it in boiling water.

STANDARDS

The **TruZone*** PFM conforms to the American Thoracic Society Standardization of Spirometry: 2005 for peak flow measurements. These standards reflect the use of 26 flow-time waveforms specifically selected to evaluate peak flow meters. Accuracy is $\pm 12\%$ or 25 L/min, whichever is greater.

① Note: The **TruZone*** PFM meets all accuracy and precision requirements that may be encountered under typical use. As with any pressure/flow rate monitor, the **TruZone*** PFM may be affected by extreme altitude variations or temperatures.** The meter is not significantly affected by relative humidity and should be used within a normal temperature range (20 to 35 degrees Celsius/68 to 95 degrees Fahrenheit). All measurements of the **TruZone*** PFM have been validated at Body Temperature and Pressure Saturated with water vapors (BTPS).

**Gardner, RM, et al. *Evaluation of Accuracy and Reproducibility of Peak Flow Meter at 1,400m*. Chest 1992; 101:948-952.

⚠ CAUTIONS

- Contact your healthcare professional immediately if you experience any unusual symptoms or if your peak flow reading indicates a worsening of your condition
- Do not leave the device unattended with children. It is not a toy.
- Device contains metal parts

📌 NOTES

- When not in use, store in a clean, dry container
- Device should be replaced after 12 months of regular use
- Follow complete cleaning instructions prior to disposal of the device
- Copies of the Instructions for Use and Daily Record sheets can be downloaded from www.truzone.ca
- Indoor use only

REORDER NUMBER

TruZone* Peak Flow Meter
(includes **TruZone*** PFM, Instructions,
ColorZone* Tapes and Daily Record) 60501070

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
PREDICTED AVERAGE PEAK EXPIRATORY FLOW (LITERS PER MINUTE)

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† Polgar G, Promadhat V. Pulmonary Function Testing in Children: Techniques and Standards. Philadelphia, W.B. Saunders Company, 1971.

السجل اليومي جـ |
DAILY RECORD OF: | REGISTRE QUOTIDIEN DE: | REGISTRO DIARIO DE: | GHI CHẾP HÀNG NGÀY VẼ:

DATE	TIME	AM	Add'l	PM
800				
750				
700				
650				
600				
550				
500				
450				
400				
350				
300				
250				
200				
150				
100				
60				

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EN Scan or contact Trudell Medical International to obtain a glossary of the symbols used on packaging and instructions for use.

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Monaghan Medical Corporation
Plattsburgh, NY 12901
USA

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