

Designed
for use
on-the-go



2 in 1 Chamber
+ Protective Case



Compact & Portable




Chambers for use at home and **on-the-go**


Always **Together™**
Inhaler +
AeroChamber® VHC

Aero2go™
Chamber


Delivers inhaler medication **better** than the inhaler alone^{1,2}



If you can **taste** your inhaler medication, it's likely not getting where it needs to go.



More medication is **delivered** to the lungs^{1,2}




Shake

inhaler immediately before use (remove cap)



Insert

inhaler into back of the chamber



Press


inhaler as you start to breathe in slowly and deeply

 **BUY NOW!**




2 in 1 Chamber + Protective Case

Your inhaler can be stored inside the chamber.




- Remove the cap from the inhaler when storing it inside the chamber
- Designed to fit most commonly prescribed inhalers
- Note: Does not fit round boot inhalers. Check compatibility before purchase.



Easy to clean + all one piece

The chamber can be used directly out of the package and should be cleaned weekly.



TRUETT MEDICAL
UK LIMITED



Designed and
Developed in Canada
with the help of
people like you

PART OF THE
AeroChamber®
FAMILY OF CHAMBERS